


Troubleshooting the Anesthesia Machine During a Procedure

By: Michelle McConnell, LVT, VTS (A&A)

Sometimes problems with the anesthesia machine are not found until the patient is under anesthesia, which is very stressful to the anesthetist. All attempts should be made to stay calm and have the ability to react quickly to solve any arising issues that could be detrimental to the patient. It is also important to rule out whether the issue is technique or machine related. There are a variety of technique issues ranging from inadequate pain control to a leaking endotracheal tube. This article will focus on equipment related issues.



A trick to testing your endotracheal tubes is to submerge the entire endotracheal tube, cuff inflated with air, into a tub of water. Verify no bubbles coming from inflation valve, pilot balloon and its line or the cuff. If bubbles are noted, discard endotracheal tube. Some silicon tubes have replacement parts available.

Reference: Page 60 Figure 7.3, Parts of the ETT, Anesthesia for Veterinary Technicians by Susan Bryant

A new leak in the anesthesia machine could present itself during a procedure. This could be discovered because the patient is not at an adequate plane of anesthesia. This happens because too much room air is entering the anesthesia machine, which dilutes the anesthetic agent. The smell of anesthetic might be present. Leaks may also create difficulty in manually ventilating the patient. Correct the leak if possible. If not, the machine should be switched out with another. If the leak is not isolated the machine should be serviced by a service company.

Another issue that may arise during a procedure is malfunctioning of the one-way valves. The purpose of the one-way valve is to ensure proper flow: Fresh gas is always flowing to the patient and expired gas is always flowing away from the patient. One-way valves can become warped or stick. If the one-way valves become warped, they will not seat correctly and will allow gas exchange. If ETCO₂ monitoring is being utilized this will reveal itself by ICO₂ levels present. Hyperventilation will occur with increased CO₂ levels and perhaps be mistaken for a light plane of anesthesia. If the one-way valves are warped new valves can be ordered and installed. If the valves are sticking, they can generally be wiped clean and put back in place. Keep breathing circuits and bags unattached to the machine when not in use to decrease the amount of condensation buildup in the anesthesia machine.



One more complication that can arise during a procedure is vaporizer output problems. If the patient is not staying at a good plane of anesthesia and is too light, the vaporizer concentration may be excessively low. Or if the patient is becoming excessively deep the vaporizer concentration may be excessively high. If either is suspected the vaporizer output should be verified by a service company or the vaporizer can be sent into a service company for cleaning and recalibration.

Works cited: Anesthesia for Veterinary Technicians by Susan Bryant