



Isoflurane Consumption Chart

Estimated CC's per Hour

Oxygen Flow (liters/min)	Vaporizer Setting					
	1.0	2.0	2.5	3.0	4.0	5.0
1.0	3.00	6.00	7.50	9.00	12.00	15.00
1.5	4.50	9.00	11.25	13.50	18.00	22.50
2.0	6.00	12.00	15.00	18.00	24.00	30.00
2.5	7.50	15.00	18.75	22.50	30.00	37.50
3.0	9.00	18.00	22.50	27.00	36.00	45.00

Sevoflurane Consumption Chart

Estimated CC's per Hour

Oxygen Flow (liters/min)	Vaporizer Setting					
	2.0	3.0	3.5	4.0	5.0	7.0
1.0	6.66	9.99	11.66	13.32	16.65	23.31
1.5	9.99	14.99	17.48	19.98	24.98	34.97
2.0	13.32	19.98	23.31	26.64	33.30	46.62
2.5	16.65	24.98	29.14	33.30	41.63	58.28
3.0	19.98	29.97	34.97	39.96	49.95	69.93

Veterinary Anesthetic
Equipment, Accessories, and Services



Visit our Web Store for easy online ordering!
www.vetamac.com
1-800-334-1583



130 Roth Court, Suite 100
 PO Box 178
 Rossville, IN 46065